Xpert Bladder Cancer Panel C104

English

INTENDED USE:

Xpert Bladder Cancer Panel C104 is intended for use as a quality control to monitor the performance of *in vitro* laboratory nucleic acid procedures for the qualitative detection of five RNA transcripts, CRH, IGF2, UPK1B, ANXA10, ABL1, performed on Cepheid GeneXpert® Systems with assays Xpert® Bladder Cancer Monitor and Xpert Bladder Cancer Detection.

The American Cancer Society estimates over 83,000 new cases of bladder cancer for 2024 and the incidence is 4 times higher in men than in women and 2 times higher in White men than in Black, Hispanic, or Asian/Pacific Islander men.¹. Timely follow-up care after treatment is extremely important for all patients because of the high likelihood of cancer recurrence or a subsequent bladder cancer.

PRODUCT SUMMARY and PRINCIPLE:

The Xpert Bladder Cancer Panel C104 kit contains Xpert Bladder Cancer Low Positive Control and Xpert Bladder Cancer Negative Control. Xpert Bladder Cancer Low Positive Control contains synthetic RNA transcripts carrying target sequences of CRH, IGF2, UPK1B, ANXA10 and ABL1 designed to be detected by Xpert Bladder Cancer assays. Xpert Bladder Cancer Negative Control contains synthetic RNA transcripts carrying target sequence of the endogenous control ABL1 only.

Routine use of quality controls that are consistent lot to lot assists the laboratory in identifying shifts, trends, and increased frequency of random errors caused by variations in the test system, such as failing reagents. Early investigation can prevent failed assay runs.

COMPOSITION:

Xpert Bladder Cancer Panel C104 contains 6 bottles of Xpert Bladder Cancer Low Positive Control and 6 bottles of Xpert Bladder Cancer Negative Control, 4mL each, composed of synthetic RNA transcripts suspended in a solution of synthetic urine, buffers, preservatives, and stabilizers. Table 1 lists the expected results for the Xpert Bladder Cancer Panel C104 when tested by Xpert Bladder Cancer assays on the GeneXpert Systems.

STORAGE and STABILITY:

Xpert Bladder Cancer Panel C104 should be stored at -25°C to -15°C. Unopened Xpert Bladder Cancer Panel C104 material is stable through the expiration date printed on the kit label when consistently stored frozen. Xpert Bladder Cancer Low Positive Control and Xpert Bladder Cancer Negative Control are for single use. Discard after use according to your local and federal regulations.

ORDERING INFORMATION

Xpert Bladder Cancer Panel C104

Part Number: C104

Kit Contains: 12 bottles x 4mL

6 each of Low Positive & Negative Controls

PRECAUTIONS and WARNINGS:

- Use Xpert Bladder Cancer Panel C104 as provided. Do not dilute or transfer to another storage tube.
- This product is intended for *in vitro* analytical testing and is provided for Research Use Only, not for use in diagnostic procedures.
- This product is clear to slightly cloudy in appearance.
- This product does not contain any biological material of human or animal origin. Universal Precautions are NOT required when handling this product.
- Reference materials should be used in accordance with local, state, federal regulations and accreditation requirements.
- Xpert Bladder Cancer Panel C104 contains guanidinium chloride which is toxic if swallowed (H302), irritating to eyes and skin (H315, H319), and is harmful to aquatic life (H402). See Regulation EC 1272/2008 (EU CLP).
- Xpert Bladder Cancer Panel C104 cannot be cloned, sold, or transferred without the explicit written consent of MMQCI.

INSTRUCTIONS FOR USE:

1. Allow the control to thaw and come to room temperature (18° – 25°C) for approximately 45 minutes. **To avoid foaming, do not use a mixer.** Control should feel warm, not cool, to the touch prior to use.
2. Use the control as provided. **DO NOT ADD TO UTR. DO NOT DILUTE.**
Note: The controls are manufactured in urine transport reagent (UTR), therefore, do not transfer the controls to UTR as this will dilute the controls and lead to false results.
3. Immediately before use, mix the control gently to avoid foaming by slowly inverting at least 10 times. Lightly tap the bottle 3 times on the bench to remove any control caught in the cap before opening the bottle.
4. Follow the assay manufacturer's Instructions for Use to process the control as a clinical sample by carefully dispensing the control directly to the sample port of the cartridge, using the transfer pipette.
5. Discard control bottle after use according to your local and federal regulations.

EXPECTED VALUES:

Xpert Bladder Cancer Panel C104 does not have assigned values. It is recommended that each laboratory establish its own performance guidelines for this product.

Table 1: Xpert Bladder Cancer Panel C104 Test Result

Kit Component:	Test Result:
Xpert Bladder Cancer Low Positive Control	Positive
Xpert Bladder Cancer Negative Control	Negative

¹ American Cancer Society. Cancer Facts & Figures 2024. Atlanta, Ga: American Cancer Society; 2024.

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