FilmArray® GI Control Panel M238 English

INTENDED USE:

The FilmArray® GI Control Panel M238 is intended for use as an external positive and negative quality control to monitor the detection and identification of *Campylobacter*, *Clostridium difficile* toxin A/B, *Plesiomonas shigelloides*, *Salmonella*, *Vibrio cholerae*, *Yersinia enterocolitica*, Enterotoxigenic *E. coli* (EPEC), Enteropathogenic *E. coli* (EPEC), Enterotoxigenic *E. coli* (ETEC) *lt/st*, Shiga-like toxin-producing *E. coli* (STEC) *stx1/stx2*, *E. coli* O157, *Shigella*/Enteroinvasive *E. coli* (EIEC), *Cryptosporidium*, *Cyclospora cayetanensis*, *Entamoeba histolytica*, *Giardia lamblia*, Adenovirus F 40/41, Astrovirus, Norovirus GI/GII, Rotavirus A, Sapovirus with the BioFire® FilmArray® Gastrointestinal (GI) Panel assay on the BioFire® FilmArray® systems. Detection of the viruses, bacteria and protozoa listed in Table 1 is an important aid to the diagnosis of gastrointestinal illness.

PRODUCT SUMMARY and PRINCIPLE:

FilmArray GI Control Panel M238 is composed of 2 controls, FilmArray GI Control M239 and FilmArray GI Control M240, both containing synthetic RNA corresponding to genome segments of pathogens listed in Table 1. The synthetic RNA is specifically designed for and intended to be used solely with the FilmArray GI Panel on the BioFire FilmArray System.

Quality controls can be used for routine monitoring of test systems, validation, verification, proficiency assessment, and training procedures. Routine use of quality controls assists the laboratory in identifying shifts, trends, and increased frequency of random errors caused by variations in the test system, such as failing reagents or malfunctioning equipment. Early investigation can prevent failed assay runs.

COMPOSITION:

The FilmArray GI Control Panel M238 is comprised of 12 single use tubes, 6 tubes of FilmArray GI Control M239 and 6 tubes of FilmArray GI Control M240, 200µL each. The controls contain synthetic RNA suspended in a non-infectious solution of buffers, preservatives, and stabilizers. Table 1 lists the pathogens that are detected in the FilmArray GI Control Panel M238 when tested by the BioFire FilmArray GI Panel.

INSTRUCTIONS FOR USE:

1. Allow the control to be tested to come to room temperature (18 °C – 25 °C), **approximately 30 minutes**.
2. Use the control as provided. **DO NOT DILUTE**.
3. Immediately before use, **mix the control thoroughly by first inverting several times followed by vortexing the tube for 3-5 seconds**. Tap the tube several times on the bench to remove any control caught in the cap before opening the tube.
4. Prepare Sample Mix according to BioFire instructions, being sure to invert the control and sample buffer mixture **at least 3 times**. Load and run a FilmArray GI Pouch, using the control as you would use a patient specimen, according to FilmArray GI Quick Guide or Instructions for Use.
5. Discard after use according to your local and federal regulations.

PRECAUTIONS and WARNINGS:

- Do not dilute. Use the control as provided.
- This product is intended for *in vitro* diagnostic use only.
- This product is only for use with BioFire FilmArray GI Panel assay on the BioFire FilmArray systems. It does not contain the entire genome of the pathogens listed in Table 1.
- Appearance: cloudy
- This product does not contain any biological material of human or animal origin. Universal Precautions are NOT required when handling this product.
- Quality control materials should be used in accordance with local, state, federal regulations and accreditation requirements.
- FilmArray GI Control Panel M238 cannot be cloned, sold, or transferred without the explicit written consent of MMQCI.

STORAGE and STABILITY:

FilmArray GI Control Panel M238 should be stored frozen (-25 °C to -15 °C). Unopened FilmArray GI Control Panel M238 is stable through the expiration date printed on the kit label when consistently stored frozen. FilmArray GI Control M239 and FilmArray GI Control M240 are for single use. Discard after use according to your local and federal regulations.

EXPECTED VALUES:

The expected results when the controls are analyzed are listed in Table 1.

Table 1: FilmArray Result Summary Interpretation

Assay	M239 Call	M240 Call
Bacteria		
<i>Campylobacter</i>	Not Detected	Detected
<i>Clostridium difficile</i> toxin A/B	Detected	Not Detected
<i>Plesiomonas shigelloides</i>	Detected	Not Detected
<i>Salmonella</i>	Not Detected	Detected
<i>Vibrio</i>	Detected	Not Detected
<i>Vibrio cholerae</i>	Detected	Not Detected
<i>Yersinia enterocolitica</i>	Detected	Not Detected
Diarrheagenic <i>E. coli</i>/Shigella		
Enterotoxigenic <i>E. coli</i> (EPEC)	Detected	Not Detected
Enteropathogenic <i>E. coli</i> (EPEC)	Ø N/A	Detected
Enterotoxigenic <i>E. coli</i> (ETEC) <i>lt/st</i>	Not Detected	Detected
Shiga-like toxin-producing <i>E. coli</i> (STEC) <i>stx1/stx2</i>	Detected	Not Detected
<i>E. coli</i> O157	Detected	Ø N/A
<i>Shigella</i> /Enteroinvasive <i>E. coli</i> (EIEC)	Detected	Not Detected
Parasites		
<i>Cryptosporidium</i>	Detected	Not Detected
<i>Cyclospora cayetanensis</i>	Not Detected	Detected
<i>Entamoeba histolytica</i>	Not Detected	Detected
<i>Giardia lamblia</i>	Not Detected	Detected
Viruses		
Adenovirus F 40/41	Detected	Not Detected
Astrovirus	Not Detected	Detected
Norovirus GI/GII	Not Detected	Detected
Rotavirus A	Not Detected	Detected
Sapovirus	Detected	Not Detected

ORDERING INFORMATION:

FilmArray® GI Control Panel M238

Part Number: M238

CONT Kit Contains: 12 tubes x 200 µL
6 each of M239 & M240

BioFire and FilmArray are trademarks of BioFire Diagnostics, LLC

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BioFire en FilmArray zijn handelsmerken van BioFire Diagnostics, LLC