## English

# Xpert BCR-ABL IS Panel C130

### INTENDED USE:

The Xpert BCR-ABL IS Panel C130 is intended for use as an assayed external quality control to monitor the performance of the *in vitro* quantitative detection of BCR-ABL1 translocation mRNA e14a2/b3a2 transcripts and the ABL1 endogenous control mRNA transcript when analyzed using the Xpert® BCR-ABL Ultra assay on Cepheid GeneXpert® Instrument Systems.

The Philadelphia chromosome, a translocation between the ABL1 gene on chromosome 9 and the BCR gene on chromosome 22, designated as t(9;22), generates the fusion gene BCR-ABL1 which is present in most chronic myelogenous leukemia patients. Quantitative monitoring of BCR-ABL1 transcripts in patient blood is an important tool for measuring response to therapy. In 2009, the World Health Organization (WHO) developed a panel of four BCR-ABL1 primary standards to establish an international scale (IS), a standardized method for reporting assay results as a ratio of fusion transcripts to control gene transcripts (%IS), useful to the harmonization of patient care across laboratories worldwide.<sup>1,2</sup> The %IS can also be expressed as molecular response (MR), the log reduction from a standardized baseline of 100% on the IS. The Xpert BCR-ABL IS Panel C130 is traceable to the WHO International Genetic Reference Panel for Quantitation of BCR-ABL Translocation (WHO Reference Panel), NIBSC code 09/138, and designed for use with the Xpert BCR-ABL Ultra assay which reports on the international scale.


### PRODUCT SUMMARY and PRINCIPLE:

Xpert BCR-ABL IS Panel C130 consists of six components. Each component contains an increasing concentration of BCR-ABL1 (e14a2/b3a2) RNA transcript mixed with a fixed concentration of ABL1 RNA transcript to produce six levels, 0.0%IS, 0.0032%IS, 0.01%IS, 0.1%IS, 1%IS and 10%IS, when analyzed on the GeneXpert system with the Xpert BCR-ABL Ultra assay. The %IS values, traceable to WHO Reference Panel, NIBSC code 09/138, are assigned to each lot of Xpert BCR-ABL IS Panel C130 according to NIBSC Instructions for Use<sup>3</sup>.

### Validation and Value Assignment

MMQCI manufactured 3 lots of Xpert BCR-ABL IS Panel C130 and tested them alongside the WHO Reference Panel, using one cartridge lot of the Xpert BCR-ABL Ultra assay. Grubb's outlier and the Bland-Altman tests were applied, lot-specific Correction Factors (CF) were calculated, and WHO-traceable %IS/MR values were assigned to each level of Xpert BCR-ABL IS Panel C130 for all 3 lots according to NIBSC Instructions for Use<sup>3</sup>. Figure 1. compares the 3 lots of Xpert BCR-ABL IS Panel C130 to the 4 members of the WHO Reference Panel. New lots of Xpert BCR-ABL IS Panel C130 will be assigned lot-specific %IS/MR values in the same manner.

**Figure 1.** Three lots of Xpert BCR-ABL IS Panel C130 calibrated to the WHO International Genetic Reference Panel for Quantitation of BCR-ABL Translocation.

### COMPOSITION:

Xpert BCR-ABL IS Panel C130 is comprised of 12 single use bottles, 2 bottles of each %IS level. The C130 bottles contain 4 mL of synthetic BCR-ABL1 RNA transcript and synthetic ABL1 control gene RNA transcript, suspended in a stabilizing matrix with a non-infectious solution of buffers and preservatives.

### STORAGE and STABILITY:

Xpert BCR-ABL IS Panel C130 should be stored at -25°C to -15°C. Unopened Xpert BCR-ABL IS Panel C130 material is stable through the expiration date printed on the kit label when consistently stored frozen. Xpert BCR-ABL IS Panel C130 components are for single use. Discard after use according to your local and federal regulations.

### INSTRUCTION FOR USE:

1. Allow the Xpert BCR-ABL IS Panel C130 component to be tested to come completely to room temperature (18°C-25°C), approximately 30 minutes, until bottles are warm to the touch.
2. Immediately before pipetting, thoroughly mix the C130 panel component by inverting 8 times followed by 2 pulse vortexes, 2-3 seconds each at maximum speed.
3. Add 4mL of the C130 panel component to 100µL of Proteinase K in a conical tube, as you would a blood specimen.
4. Continue with the assay procedure according to manufacturer's instructions.
5. Discard after use according to local and federal regulations.

### PRECAUTIONS and WARNINGS:

- Use the control as provided. Do not dilute or transfer to another storage tube.
- This product is intended for *in vitro* diagnostic use only.
- Use Xpert BCR-ABL IS Panel C130 only with Xpert BCR-ABL Ultra assay.
- Xpert BCR-ABL IS Panel C130 is not intended to be used for calibration of the Xpert® BCR-ABL Ultra assay.
- This product is slightly cloudy in appearance.
- This product does not contain any biological material of human or animal origin. Universal Precautions are NOT required when handling this product.
- Xpert BCR-ABL IS Panel C130 cannot be cloned, sold, or transferred without the explicit written consent of MMQCI.

### EXPECTED VALUES:

Locate the appropriate WHO-traceable %IS/MR values assigned to your lot of Xpert BCR-ABL IS Panel C130 on the Data Sheet found in each kit box of Xpert BCR-ABL IS Panel C130. It is important to notice that the WHO-traceable values were assigned by testing with one lot of Xpert BCR-ABL Ultra cartridges. Each laboratory should establish their own %IS/MR ranges. Linearity can be confirmed by performing a linear regression with an expected correlation coefficient (R<sup>2</sup>) at or above 0.9. Please refer to the Xpert BCR-ABL Ultra Package Insert for expected assay performance specifications.

Routine use of quality controls that are consistent lot-to-lot assists the laboratory in identifying shifts, trends, and increased frequency of random errors caused by variations in the test system.

### ORDERING INFORMATION:

Xpert BCR-ABL IS Panel C130

**Part Number: C130**

**CONT** Kit Contains: 12 bottles x 4mL; 2 x 0%IS, 2 x 0.0032%IS, 2 x 0.01%IS, 2 x 0.1%IS, 2 x 1%IS, 2 x 10%IS

### References

- <sup>1</sup>Branford S et al. Desirable performance characteristics for BCR-ABL measurement on an international reporting scale to allow consistent interpretation of individual patient response and comparison of response rates between clinical trials. *Blood* 2008, 112:3330-38
- <sup>2</sup>White HE et al. Establishment of the first World Health Organization International Genetic Reference Panel for quantitation of BCR-ABL mRNA. *Blood* 2010, 116:e111-117
- <sup>3</sup>1st WHO International Genetic Reference Panel for Quantitation of BCR-ABL Translocation, NIBSC code: 09/138 Instructions for use (Version 4.0, Dated 13/12/201)

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### Referenties

- Branford S et al. Desirable performance characteristics for BCR-ABL measurement on an international reporting scale to allow consistent interpretation of individual patient response and comparison of response rates between clinical trials. *Blood* 2008, 112:3330-38
- White HE et al. Establishment of the first World Health Organization International Genetic Reference Panel for quantitation of BCR-ABL mRNA. *Blood* 2010, 116:e111-117
- 1st WHO International Genetic Reference Panel for Quantitation of BCR-ABL Translocation, NIBSC code: 09/138 Instructions for use (versie 4.0, gedateerd 13/12/2012)