



## SPOTFIRE® RSP Positive Control

### INTENDED USE:

SPOTFIRE® RSP Positive Control is intended for use (as applicable) as a reference material to monitor performance of *in vitro* laboratory nucleic acid testing procedures for the qualitative detection of *Bordetella parapertussis*, *Bordetella pertussis*, *Chlamydia pneumoniae*, *Mycoplasma pneumoniae*, *Streptococcus dysgalactiae* (Group C/G Strep), *Streptococcus pyogenes* (Group A Strep), Adenovirus, seasonal Coronavirus, Coronavirus SARS-CoV-2, Human Metapneumovirus, Human Rhinovirus/Enterovirus, Influenza A Subtype H1-2009, Influenza A Subtype H3, Influenza B, Parainfluenza Virus and Respiratory Syncytial Virus using the BIOFIRE® SPOTFIRE Respiratory (R) Panel (SPOTFIRE R Panel), BIOFIRE SPOTFIRE Respiratory (R) Panel Mini (SPOTFIRE R Panel Mini) and BIOFIRE SPOTFIRE Respiratory Sore Throat Panel (SPOTFIRE R/ST Panel) assays performed on the BIOFIRE SPOTFIRE System. SPOTFIRE RSP Positive Control is comprised of *in vitro* RNA transcripts and stabilizing solution and is designed for and intended to be used solely with the SPOTFIRE R Panel, SPOTFIRE R Panel Mini and SPOTFIRE R/ST Panel assays. This product is not intended to replace manufacturer internal controls provided with the device.

Quality control materials should be used in accordance with local, state, federal regulations, and accreditation requirements.

### PRODUCT SUMMARY and PRINCIPLE:

SPOTFIRE RSP Positive Control contains surrogate control material composed of synthetic RNA corresponding to genome segments of pathogens listed as detected in Tables 1, 2 and 3.

Quality controls can be used for routine monitoring of test systems, validation, verification, proficiency assessment, and training procedures. Routine use of quality controls assists the laboratory in identifying shifts, trends, and increased frequency of random errors caused by variations in the test system, such as failing reagents or malfunctioning equipment. Early investigation can prevent failed assay runs.

### COMPOSITION:

SPOTFIRE RSP Positive Control kit is comprised of 6 single use tubes, 300µL each, containing synthetic RNA suspended in a non-infectious solution of buffers, preservatives and stabilizers. The synthetic RNA is specifically designed for and intended to be used solely with the SPOTFIRE R Panel, SPOTFIRE R Panel Mini and SPOTFIRE R/ST Panel assays on the BIOFIRE SPOTFIRE System. Tables 1 and 2 list the pathogens that are detected in the SPOTFIRE RSP Positive Control when tested with the SPOTFIRE R Panel, SPOTFIRE R Panel Mini, and Table 3 for SPOTFIRE R/ST Panel assays respectively.

### INSTRUCTIONS FOR USE:

- The controls may be shipped on dry ice or cold packs and placed at room temperature (18° – 25°C), or refrigerated (2° – 8°C) upon receipt.
- Allow the control to come completely to room temperature (18° – 25°C) before use, **approximately 30 minutes** if frozen due to dry ice shipping.
- Use the control as provided. Do Not Dilute.
- Immediately before use, mix the control **thoroughly** by first flicking the tubes 3 times followed by inverting at least 3 times. Tap the tube 3 times on the bench to remove any control caught in the cap before opening the tube.
- Refer to the Quality Control Testing section of the SPOTFIRE R Panel, SPOTFIRE R Panel Mini, or SPOTFIRE R/ST Panel Quick Guide to prepare, hydrate and load the SPOTFIRE R Panel Pouch, SPOTFIRE R Panel Mini Pouch, or SPOTFIRE R/ST Panel Pouch and to test the QC sample.  
*Note: It is important to **thoroughly** mix the QC sample with the Sample Buffer by inverting the Sample Injection Vial 3 times, as described in the Quick Guide.*
- Discard after use according to your local and federal regulations.

For complete SPOTFIRE R Panel, SPOTFIRE R Panel Mini or SPOTFIRE R/ST Panel instructions, refer to the appropriate Quick Guide or Instruction Booklet provided by BioFire Diagnostics, LLC.

### STORAGE and STABILITY:

SPOTFIRE RSP Positive Control can be shipped on dry ice or cold packs, and stored at room temperature (18° – 25°C) or refrigerated (2° – 8°C) upon receipt. The controls are stable through the expiration date printed on the kit label. SPOTFIRE RSP Positive Control kit components are for single use. Discard after use according to your local and federal regulations.

### PRECAUTIONS and WARNINGS:

- Do not dilute. Use the control as provided.
- This product is intended for *in vitro* analytical testing and is provided for Research Use only, not for use in diagnostic procedures.
- SPOTFIRE RSP Positive Control is for use with the SPOTFIRE R Panel, SPOTFIRE R Panel Mini or SPOTFIRE R/ST Panel assays on the BIOFIRE SPOTFIRE System. It does not contain the entire genome of pathogens listed in Table 1.
- Appearance: Positive control is slightly cloudy.
- This product does not contain any biological material of human or animal origin. Universal Precautions are NOT required when handling this product.
- SPOTFIRE RSP Positive Control cannot be cloned, sold, or transferred without the explicit written consent of MMQCI.

### EXPECTED VALUES:

The expected results when the SPOTFIRE RSP Positive Control is tested as a Quality Control (QC icon) with SPOTFIRE R Panel (Table 1), SPOTFIRE R Panel Mini (Table 2) and SPOTFIRE R/ST Panel (Table 3).

**Table 1. SPOTFIRE R Panel Results**

<span style="border: 1px solid black; border-radius: 50%; padding: 2px;">R</span> Positive QC Test	
<b>PASS</b>	
<b>Bacteria</b> <i>Positive ✓ Bordetella parapertussis</i> <i>Positive ✓ Bordetella pertussis</i> <i>Positive ✓ Chlamydia pneumoniae</i> <i>Positive ✓ Mycoplasma pneumoniae</i>	<b>Viruses</b> <i>Positive ✓ Adenovirus</i> <i>Positive ✓ Coronavirus SARS-CoV-2</i> <i>Positive ✓ Coronavirus (seasonal)</i> <i>Positive ✓ Human Metapneumovirus</i> <i>Positive ✓ Human Rhinovirus/Enterovirus</i> <i>Positive ✓ Influenza A Virus</i> <i>Positive ✓ Influenza A virus A/H1-2009</i> <i>Positive ✓ Influenza A virus A/H3</i> <i>Positive ✓ Influenza B Virus</i> <i>Positive ✓ Parainfluenza Virus</i> <i>Positive ✓ Respiratory Syncytial Virus</i>

**Table 2. SPOTFIRE R Panel Mini Results**

<span style="border: 1px solid black; border-radius: 50%; padding: 2px;">R</span> Positive QC Test	
<b>PASS</b>	
<b>Viruses</b> <i>Positive ✓ Coronavirus SARS-CoV-2</i> <i>Positive ✓ Human Rhinovirus</i> <i>Positive ✓ Influenza A Virus</i> <i>Positive ✓ Influenza B Virus</i> <i>Positive ✓ Respiratory Syncytial Virus</i>	

**Table 3. SPOTFIRE R/ST Panel Results**

<span style="border: 1px solid black; border-radius: 50%; padding: 2px;">R</span> <span style="border: 1px solid black; border-radius: 50%; padding: 2px;">ST</span> Positive QC Test	
<b>PASS</b>	
<b>Bacteria</b> <i>Positive ✓ Bordetella parapertussis</i> <i>Positive ✓ Bordetella pertussis</i> <i>Positive ✓ Chlamydia pneumoniae</i> <i>Positive ✓ Mycoplasma pneumoniae</i> <i>Positive ✓ Streptococcus dysgalactiae (group C/G Strep)</i> <i>Positive ✓ Streptococcus pyogenes (group A Strep)</i>	<b>Viruses</b> <i>Positive ✓ Adenovirus</i> <i>Positive ✓ Coronavirus SARS-CoV-2</i> <i>Positive ✓ Coronavirus (seasonal)</i> <i>Positive ✓ Human Metapneumovirus</i> <i>Positive ✓ Human Rhinovirus/Enterovirus</i> <i>Positive ✓ Influenza A Virus</i> <i>Positive ✓ Influenza A virus A/H1-2009</i> <i>Positive ✓ Influenza A virus A/H3</i> <i>Positive ✓ Influenza B Virus</i> <i>Positive ✓ Parainfluenza Virus</i> <i>Positive ✓ Respiratory Syncytial Virus</i>

BIOFIRE and SPOTFIRE are trademarks of BioFire Diagnostics, LLC

### ORDERING INFORMATION:

Product Name: SPOTFIRE RSP Positive Control  
 Part Number: M42638  
 Kit Contains: 6 tubes of Positive control x 300µL