

Development of a *CFTR* Control Panel to Monitor All ACMG-100 Variants

Jessica Hayden¹, Aldona Karaczyn¹, Joan Gordon¹, **Joe Salem¹**, Jacob Schurkman¹, Tania Spenlinhauer¹, Andrew Bradford², Katarzyna Gawlik², Laura Mendez²

¹Maine Molecular Quality Controls Inc., Saco, Maine

²Agena Bioscience, San Diego, California


G021

Introduction

In 2001, the American College of Medical Genetics and Genomics (ACMG) published a list of 23 variants that should be included in carrier screening for cystic fibrosis. In 2023, ACMG updated their recommended variant list to 100 [1]. Next-generation sequencing (NGS) and MassARRAY are now routine methods to detect *CFTR* variants in clinical laboratories. Although testing is routine, template preparation and bioinformatics pipelines are complex. Sequence noise and amplification errors must be reliably distinguished from true variant calls. Sanger sequenced-verified external controls are critical to monitor the ability of any NGS or MassARRAY test system to correctly identify both simple and complex genome variants. Here, we describe development and verification of a safe, non-infectious, and easy-to-use 6-bottle *CFTR* control panel containing all 100 ACMG-recommended variants, as well as variants challenging to NGS systems.

Materials and Methods

Twenty-eight additional variants were added to existing plasmids from the previous version of the control panel (see Table 2, asterisk*). Plasmids contain all *CFTR* exons with intronic borders and were used to create a synthetic panel containing a total of 211 *CFTR* variants and all 100 of the recommended ACMG variants. Fragments containing *CFTR* variants were designed *in-silico*, and incorporated using Gibson assembly. Assembled constructs were transformed into bacteria and screened using Oxford Nanopore sequencing. Plasmids were purified from the selected clones and bi-directionally sequenced using Sanger (Figure 1).



Bottle	Control Name
A	G211Aplus
B	G211Bv1.1
C	G211Cplus
D	G211Dv1.1
E	G211Eplus
F	G211Fplus

Table 1 (left). NGS CF Control G211plus panel configuration:

All 211 *CFTR* variants are split across 6 controls. Purified plasmids each containing unique *CFTR* variants were mixed at a 1:1 ratio to generate heterozygous calls.

Three unique lots (6 control bottles per lot) were tested by NGS using the Illumina TruSight™ CF Clinical Sequencing Assay on the MiSeqDx at MMQCI. One lot of G211plus controls was also provided to Agena Bioscience for testing on the new *CFTR*100+ assay using the MassARRAY® system. *CFTR*100+ is designed to detect all ACMG-100 variants. Figure 2 (below) shows a brief description of their workflows.




Table 2: List of all 211 variants in the NGS CF Control G211plus panel (ACMG-100 is **bolded**, the 28 newly added ACMG-100 variants have an asterisk*)

	cHGVs (Legacy)		
c.1A>G (M1V)	c.805_806delAT (936delTA)	c.2260G>C (V754M)	c.3454G>C (D1152H)
c.4C>T (Q29X)*	c.868C>T (Q290X)	c.290C>T (R764K>C>T)	c.3472C>T (R1158X)
c.91C>T (R31C)	c.929TCT[2] (delF311)	c.1584+1G>A (T17+1G>A)*	c.3484C>T (R1162X)
c.115C>T (Q39X)	c.948delT (I078delT)	c.1585+1G>A (T17+1G>A)*	c.3485G>T (R1162L)
c.164T>A (D36>2T>A)	c.980C>T (G338X)	c.1585delG (T17+1G>G>A)	c.3520delC (S609delC)
c.178G>C (E60X)	c.1000C>T (R334V)	c.1624G>T (G542X)	c.3535_3536delCCAA (S667del4)
c.200C>T (P67L)	c.1007T>A (I336K)	c.1645A>C (S549R)	c.3587C>G (S1196K)
c.232C>T (R75X)	c.1013C>T (T338I)	c.1646G>A (S549N)	c.3611G>A (W1204X (S743G>A))
c.254G>A (G58R)	c.1021T>C (S341P)	c.1647T>C (S549K)	c.3612G>A (W1204X (CF746G>A))
c.262_263delTT (294delTT)	c.1021_1022dupCT (T154insTC)	c.1651G>A (G551S)*	c.3659delC (S791delC)
c.271G>A (G91R)*	c.1029delT (I161delC)*	c.1652G>A (G551D)	c.3700A>G (I1234V)
c.273-1G>A (A05+1 G>A)	c.1040G>A (R347H)	c.1654C>T (Q552X>T)	c.3705T>G (S1233R)
c.273+3A>C (A05+3A>C)	c.1040C>C (R347P)	c.1657+2_2657+3insA (2789+2insA)	c.3717+5G>A (S849+5G>A)*
	c.1055G>A (R352D)	c.1672T>C (L558S)	c.3718-247C>T (3849+108k C>T)
	c.1081delT (I213delT)	c.1675G>A (A559T)	c.3717+4A>G (S849+4A>G)
c.274+1G>A (delC+G>A)	c.1090T>C (S364P)	c.1679C>C (R500T)	c.3737C_2738insG (2669insG)
c.274G>A (P191Y)*	c.1116+1G>A (L248+1G>A)	c.1679G>A (R500K)	c.3744delA (S761delA)
c.274G>T (E92X>G>T)	c.1130delA (I259insA)	c.1679+1G>A (R11+1G>A)*	c.3752G>A (S1251N)
c.292C>T (Q98X)	c.1155_1156dupTA (I288insTA)	c.1680-886A>G (1811+1.6 kb)	c.3807A>G (S1255X>A>G)
c.293A>G (Q98R)*	c.1202G>A (W401X (TAG))	c.1680A>C (R505R)*	c.3764C>A (S1255X>A)
c.313delA (444delA)	c.1206G>A (W401X (TGA))	c.1680C>C (Q970R)	c.3773delT (S905delT)
c.325_327delinsG (457TAT>G)	c.1209+1G>A (L341+1G>A)	c.1682C>A (G970M)	c.3808del (294delG)*
c.328G>C (D110H)	c.1327_1330dupGATA (I461ins4)	c.2930C>T (S977L)	c.3808G>A (D1270N)
c.349C>T (R117C)	c.1344C>A (A455L)	c.2980G>A (S128G>A)	c.3846G>A (W1282X)
c.358G>A (R117H)	c.1367T>C (V45A)*	c.3272C>G (G575A)	c.3873C>G (A05+1G>A)
c.366T>A (Y122X)	c.1373del (L504delG)*	c.1753G>T (E585X)	c.3889dup (4016insT)
c.442delA (L574delA)	c.1393-1G>A (I525+1G>A)	c.2991G>C (L997E (G>C))	c.3900C>G (N1303K)
c.443T>C (I48T)	c.1397C>A (S460X (C>A))	c.3039delC (T11delC)	c.3917C>T (Q1313X)
c.489+1G>T (G21+1 G>T)	c.1397C_3046X (C>G)	c.1397C_3046X (C>G) (CTG919del4)	c.4077_408delA (A4209GTT>AA)
	c.1401T>C (L467P)	c.1837C>C (I1027T)	c.4251delA (G382delA)
c.532G>A (G178R)	c.1408G>A (M470V)	c.1841A>G (D614G)	c.54-5940_273+1025del1 (CFTRdele2,3)
c.571T>G (P191Y)*	c.1418delA (S184delA)	c.1695G>A (G520K)	c.1418-26A>G (272+26A>G)
c.579+1G>T (T11+1 G>T)	c.1438G>T (G480C)	c.1882C>A (G628R)	c.1947C>T (L1065P)
c.579+3A>G (T11+3 A>G)	c.1466C>A (S489X)	c.1923_1931delinsA (2055del89>A)	c.3196C>T (R1066G)
c.579+5G>A (T11+5 G>A)	c.1475C>T (S492I)	c.1976delA (I208delA)	c.3197G>A (R1066I)
c.580+1G>T (T12+1 G>T)	c.1477C>T (Q493X)	c.2002C>T (R666C)	c.3206C>T (R1070N)
c.595C>T (I1199Y)	c.1519A>G (S506V)	c.2012delT (I2143delT)	c.3209G>A (R1070Q)
c.613C>T (P205S C>T)	c.1519A>G (S507V)	c.2051_2052delinsG (2183AA>G)	c.3230T>C (L1077P)
c.617T>G (L286W)	c.1516ATC[1] (S507del)	c.2052dupA (I2184insA)	c.3266G>A (W1089X)
c.637A>C (L181X)*	c.1521_1523delCTT (F508del)	c.2053delA (I2184delA)	c.3276C>A (Y1092X (C>A))
c.658C>T (Q220X)	c.1523T>G (F508K)	c.2125C>T (R709X>C>T)	c.3276G>G (Y1092X (C>G))
c.680T>C (L227R)	c.1545_1546delTA (I677delTA)	c.2128A>T (R710X)	c.3294G>C (W1098C*)
c.695T>A (V232D)	c.1558G>T (V520R)	c.2175dupA (2307insA)	c.3302T>A (M1101K)
c.723_743-1del (853del22)	c.1572X>A (C534X)*	c.3310C>T (R1104X)	c.3310T>G (L1732X)
c.803delA (935delA)		c.2215delAG (L241delG)	c.3355C>T (R1118P*)

Results (Next Generation Sequencing)

Bottle	P/N	Lot #	Total (n) across 3 sequencing runs	NGS Call rate	Correct calls/Expected Calls	% Concordance
A	NGS CF Control Panel G211Aplus	A27MAY25A	6	100.00%	38/38	100.0%
		E18JUN25A	6	100.00%	38/38	100.0%
		E02JUL25A	6	100.00%	38/38	100.0%
B	NGS CF Control Panel G211Bv1.1	B28MAY25A	6	99.04%	38/38	100.0%
		C23JUN25A	6	99.04%	38/38	100.0%
		B08JUL25A	6	99.04%	38/38	100.0%
C	NGS CF Control Panel G211Cplus	C27MAY25A	6	99.65%	36/36	100.0%
		G18JUN25A	6	99.65%	36/36	100.0%
		H02JUL25A	6	99.65%	36/36	100.0%
D	NGS CF Control Panel G211Dv1.1	D28MAY25A	6	99.04%	39/39	100.0%
		E23JUN25A	6	99.04%	39/39	100.0%
		D08JUL25A	6	99.04%	39/39	100.0%
E	NGS CF Control Panel G211Eplus	E28MAY25A	6	99.94%	33/33	100.0%
		G23JUN25A	6	99.94%	33/33	100.0%
		E08JUL25A	6	99.94%	33/33	100.0%
F	NGS CF Control Panel G211Fplus	F28MAY25A	6	99.04%	44/45	99.8%
		J23JUN25A	6	99.04%	44/45	99.8%
		U09JUL25A	6	99.04%	44/45	99.8%
Total n = 108				Avg call rate = 99.45%	684/687	99.60%

P/N	G211Aplus	G211Aplus	G211Aplus	G211Bv1.1	G211Bv1.1	G211Bv1.1
Lot	A27MAY25A	E18JUN25A	E02JUL25A	B28MAY25A	C23JUN25A	B08JUL25A
N replicates	6	6	6	6	6	6
N calls	216	216	216	222	222	222
Mean AF	0.499	0.509	0.508	0.477	0.515	0.469
Median AF	0.498	0.518	0.516	0.468	0.505	0.455
Minimum AF	0.316	0.359	0.346	0.385	0.442	0.405

P/N	G211Cplus	G211Cplus	G211Cplus	G211Dv1.1	G211Dv1.1	G211Dv1.1
Lot	C27MAY25A	G18JUN25A	H02JUL25A	D28MAY25A	E23JUN25A	D08JUL25A
N replicates	6	6	6	6	6	6
N calls	210	210	210	228	228	228
Mean AF	0.506	0.509	0.509	0.488	0.484	0.450
Median	0.487	0.497	0.492	0.479	0.471	0.435
Minimum	0.366	0.396	0.401	0.347	0.378	0.340

P/N	G211Eplus	G211Eplus	G211Eplus	G211Fplus	G211Fplus	G211Fplus
Lot	E28MAY25A	G23JUN25A	E08JUL25A	F28MAY25A	J23JUN25A	U09JUL25A
N replicates	6	6	6	6	6	6
N calls	192	192	192	252	252	252
Mean AF	0.503	0.495	0.430	0.473	0.457	0.475
Median	0.509	0.497	0.429	0.472	0.456	0.476
Minimum	0.338	0.357	0.304	0.338	0.342	0.359

Table 3: Breakdown of variant types and subtypes in the control panel.

Variant Types & Subtypes	Number of variant type in the control panel
Deletions	30
1 bp	23
2-6 bp	6
22 bp	1
Deletion-Insertion (Delins)	4
Duplications & Insertions	10
1 bp	7
2-4 bp	3
Large Genomic Rearrangements	2
Tandem repeat variants	4
Substitutions (exon)	131
A>C	2
A>G	7
A>T	1
C>A	8
C>G	4
C>T	33
G>A	33
G>C	8
G>T	11
T>A	6
T>C	10
T>G	8
Substitution (intronic)	30
A>C	1
A>G	5
C>T	1
C>A	18
G>T	4
T>A	1

Table 4. Total NGS sample replicates and percent correct calls: A total of 18 G211plus lots (6 controls x 3 lots each) were manufactured and tested across 3 Illumina TruSight CF runs. The Clinical Sequencing Assay reports detected variants along with their specific allele frequencies. Call rate, a reported percentage, is the number of variant regions that meet a defined confidence value threshold divided by the total regions targeted. Each control lot was tested 6 times (total n = 108). Bottles A through E had 100% correct calls across all NGS runs. Variant c.3700A>G in Bottle F is not reported because its allele frequency is below a proprietary cutoff set by Illumina. The low frequency is caused by a known near-neighbor interference issue with c.3612G>A. c.3700A>G can still be observed using IGV at ~15% allele frequency (Figure 3, below).

Figure 3. IGV Visualization of c.3700A>G in G211Fplus. The left image is sorted by base at c.3700 and shows the frequency of each call. The right image highlights c.3700A>G in paired, unsorted reads.




Table 5. Lot-to-lot reproducibility of allele frequency: Allele frequency (AF) summary statistics, including averages, medians, and minimums in each control were calculated using JMP 18. Average allele frequencies ranged from 0.430 to 0.515. Note that PolyT/PolyTG variants and the two LGRs (CFTRdele2,3 and CFTRdele22,23) are reported without frequencies. Homozygous variant c.1408G>A is present in Bottles B and D.

Table 6. Examples with IGV images highlighting the importance of quality controls when interpreting bioinformatic results. Below are examples of how bioinformatics pipelines and variant location can impact NGS reporting.

Variant	Control	TruSight AF Skew	Cause of Skew	Details
c.3302T>A	G211Aplus	Low	Illumina bioinformatics	The presence of the T>A substitution, which is adjacent to an 'A'-homopolymer causes a misalignment at the end of reads. The shift is counted as a "deletion" rather than the 'A' allele.
c.1013C>T	G211Eplus	Low	Illumina bioinformatics	When present as the last base at the end of reads, the allele is removed by bioinformatic soft-trimming which significantly reduces the frequency.
c.1516A>G	G211Cplus	High	Overlapping variants in the control	c.1516A>G reports at near homozygous frequency because it is paired with deletion c.1516ATC[1] in the same control. One allele is 'G' and the other allele is deleted.
c.1021_1022dupTC	G211Aplus	High	Illumina bioinformatics	Allele frequency is elevated resulting from a population of unpaired reads containing the TC-duplication.





Figure 4. Results from G211plus testing on the Agena CFTR100+ assay: (A) Variant calling heat map of G211plus controls and Coriell samples. (B) Allelic ratios across a dilution series of the control material. The controls were originally developed for NGS. The dilution series suggests that the *CFTR*100+ assay is more sensitive than NGS. (C) Box plot of assay yield across a dilution series. Yield is a metric that indicates assay performance. Yield > 0.5 is passing. Yield > 0.8 is the target performance goal. G211plus controls generated a median yield of 0.98. (D) Causes for G211plus discrepant calls when tested on the *CFTR*100+ assay.

Conclusions

- We have developed a new *CFTR* external control panel that can be used to monitor all 100 of the updated ACMG-recommend variants for carrier screening.
- The control panel includes clinically relevant variants of many types that challenge bioinformatic and variant reporting pipelines, including 161 substitutions representing all possible nucleotide changes, 30 small deletions, 4 complex insertion/deletion variants, 10 small insertions/duplications, 2 large genomic rearrangements, and Poly-T/Poly-TG repeats.
- The control panel serves to confidently control for NGS and MassARRAY wet-lab methods despite their fundamentally different chemistries, detection, and output.
- G211plus controls ensure bioinformatic pipelines can detect all relevant *CFTR* variants and reveal causes of allelic skew from bioinformatic processing.

Acknowledgements: Special thanks to Agena Bioscience and their colleagues: Andy Bradford for the MassARRAY figures, Laura Mendez who ran the titration experiment, and Katarzyna Gawlik who tested the Coriell samples.

[1] Deignan JL, Gregg AR, Monaghan KG, et al. Updated recommendations for CFTR carrier screening: A position statement of the American College of Medical Genetics and Genomics. *Genetics in Medicine*. 2023