TECHNICAL NOTE # 17-001

RE: NGS CF Control Panel G211 v1.1 (p/n: G211v1.1)
ASSAY: TruSight™ Cystic Fibrosis 139-Variant Assay (Illumina)
DATE: 10/25/2022

This Technical Note serves as a guide when using NGS CF Control Panel G211 v1.1 with TruSight Cystic Fibrosis 139-Variant Assay. The protocol below is to be used along with the Instructions for Use in the NGS CF Control Panel G211 v1.1 Package Insert and the TruSight Cystic Fibrosis 139-Variant Assay Package Insert. Table I lists the genotype calls expected to be seen in the TruSight Cystic Fibrosis 139-Variant Assay report.

NGS CF Control Panel G211 v1.1 Protocol for use with the TruSight Cystic Fibrosis 139-Variant Assay:

Controls are ready to use; no extraction is needed.

1. Thaw all controls completely.
2. Allow the controls to be tested to come to room temperature (18° – 25°C).
3. Mix well by flicking the control tubes and vortex immediately before using.
4. Quick spin to remove droplets from the tube cap before opening.
5. Add 5 µL of each control to the appropriate wells in the **HYB** plate of Illumina’s TruSight Cystic Fibrosis 139-Variant Assay (as described in the TruSight Cystic Fibrosis 139-Variant Assay protocol). The controls cannot be quantified, they are manufactured to have approximately equivalent copy numbers of the target gene as genomic DNA.
6. Follow all other procedures as outlined in TruSight Cystic Fibrosis 139-Variant Assay protocol.
7. Opened material should be tightly capped and returned to freezer ($\leq -20^{\circ}\text{C}$) immediately after use.

Table I on the following page lists the expected calls for the NGS CF Control Panel G211 v1.1 when tested with the Illumina TruSight Cystic Fibrosis 139-Variant Assay. Please consult with Maine Molecular Quality Controls, Inc. if assistance is required to get the correct genotypes listed in the table below.

Table I Footnotes:

1. PolyTG/PolyT variants are only reported if the R117H variant is present.
2. Although the three benign polymorphisms, I506V, I507V and F508C, are present within the NGS CF Control Panel G211 v1.1, these variants will not be reported because the TruSight Cystic Fibrosis 139-Variant Assay reports them “conditionally,” only when F508del or I507del is homozygous mutant.

General Notes:

- Expected call rate is approximately 100%. If Call Rate falls below 99%, it is recommended to contact MMQCI for further guidance.
- Samples **cannot** be quantified or measured by spectrophotometry or fluorometry. Synthetic DNA does not have enough mass to be quantified by these methods; however, the controls are manufactured to contain CFTR target sequence at copy numbers approximately equivalent to that of genomic DNA.
- All 139 variants detected by the TruSight Cystic Fibrosis 139-Variant Assay are present in the NGS CF Control Panel G211 v1.1.

Please call us at 207-885-1072 x 1 if you have any questions.

Sincerely,
Maine Molecular Quality Controls, Inc.
www.mmqci.com

TECHNICAL NOTE # 17-001

Table I. NGS CF Control Panel G211 v1.1 Expected Calls with the TruSight Cystic Fibrosis 139-Variant Assay

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Bold = ACMG-23
Italics = Conditionally Reported
N/A = Mutation not present
N/R = Not Reported