



Quality Controls for Infectious Disease Testing on the Cepheid GeneXpert[®]



DEPENDABLE

Third party independent controls with lot to lot consistency tested over multiple sites, instruments, operators, and cartridge lots



NON-INFECTIOUS

Do not contain any biological material of human or animal origin. Precautions are not required when handling the products



EFFICIENT

Full process DNA and/or RNA controls that mimic patient sample to monitor analytical performance of the extraction, amplification and detection steps

Routine use of quality controls can assist in identifying shifts, trends, and increased frequency of random errors caused by variations in the test system, such as failing reagents or malfunctioning equipment. Early investigation can prevent failed assay runs.





Tuberculosis

INTROL® TB Panel M114

Catalog No.: M114-5

Kit Contains: 5 bottles each of TBWT-04, TBMDR1-04, & TBNEG x 1 mL

Catalog No.: M114-10

Kit Contains: 10 bottles each of TBWT-04, TBMDR1-04, & TBNEG x 1 mL

INTROL® TB Panel M114plus

Catalog No.: M114plus

Kit Contains: 3 bottles each of TBWT-04, TBMDR1-04, TB1081MDR2, & TBNEG x 1 mL

INTROL® TBNEG-30

Catalog No.: M115

Kit Contains: 30 bottles of TBNEG x 1 mL

INTROL® TB-XDR Panel M345

Catalog No.: M345-5

Kit contains: 5 bottle each of TB-XDR WT, TB-XDR MUT, & TB-XDR NEG x 1 mL



Respiratory

Xpert SARS-CoV-2 Controls Panel M430

Catalog No.: M430

Kit Contains: 3 tubes each of Positive & Negative x 300 µL

Birlinn® CoV2-Flu-RSV Control Panel M509*

Catalog No.: M509

Kit Contains: 6 tubes each of Positive & Negative x 800 µL

**Tested and compatible with Xpert® Xpress Flu, Xpert® Xpress CoV-2 plus, and Xpert® Xpress CoV-2/Flu/RSV plus assays*

NEW



Antibiotic Resistance

Xpert Carba-R QC Panel M219

Catalog No.: M219

Kit Contains: 6 tubes each Positive & Negative x 50 µL

Maine Molecular Quality Controls, Inc. | 23 Mill Brook Road, Saco, ME 04072, USA

www.mmqci.com

info@mmqci.com

orders@mmqci.com

207-885-1072

207-885-1079