



# Infectious Disease Testing on the Diasorin LIAISON PLEX®



## MULTIPLEX

Monitors the amplification and detection of all 19 pathogens included in the LIAISON PLEX® Respiratory Flex Assay with only 2 tubes, saving costs of cartridges used.



## READY & SAFE TO USE

No reconstitution or pooling required. Contains carefully crafted genome segments that are 100% non-infectious and safe for the lab to handle.



## TRUE NEGATIVE

Includes a true negative control that contains no target template, able to detect low-levels of contamination.



## VERSATILE

Can be used for validation, verification, routine monitoring, and operator training, easy to include in your IQCP.



## Respiratory

### LIAISON PLEX RP Control Panel M518

Catalog No.: M518

Kit Contains: 12 tubes x 300 µL; 4 each of Pos A, Pos B, & Negative controls



## Blood Culture

### LIAISON PLEX Yeast Control Panel

Catalog No.: M544

Product in Development

### LIAISON PLEX GN Control Panel

Catalog No.: M548

Product in Development

### LIAISON PLEX GP Control Panel

Catalog No.: M552

Product in Development



Maine Molecular Quality Controls, Inc. | 23 Mill Brook Road, Saco, ME 04072, USA

 [www.mmqci.com](http://www.mmqci.com)

 [info@mmqci.com](mailto:info@mmqci.com)

 [orders@mmqci.com](mailto:orders@mmqci.com)

 207-885-1072

 207-885-1079

MBR076-03.26.2026