

Reference Materials for Oncology Testing on the Cepheid GeneXpert® System



Ensure Reliable Cancer Diagnostics

Accurate detection and monitoring is crucial for early diagnosis, prognosis, and treatment of cancer patients. MMQCI's Xpert multiplex quality controls are uniquely designed to be carried through all phases of molecular testing on the Cepheid GeneXpert Systems.

DEPENDABLE

Third party independent controls tested over multiple lots, sites, operators, and instruments

SAFE

Do not contain any biological material of human or animal origin. Universal precautions are not required when handling the products

PREVENTATIVE

Help to identify shifts, trends, and increased frequency of random errors caused by variations in the test system

READY TO USE

Convenient single use 4mL format to be used as provided, no dilution, reconstitution, or pooling needed

LONG SHELF LIFE

Can be stored between -25°C and -15°C for up to 36 months from date of manufacture

VERSATILE

Can be used for routine monitoring, validation, verification, proficiency assessment, range confirmation and training



Bladder

Xpert Bladder Cancer Panel C104

Catalog No.: C104

Kit Contains: 6 bottles each of Low Positive & Negative controls x 4 mL



BCR-ABL

Xpert BCR-ABL IS Panel C130*

Catalog No.: C130

Kit contains: 2 bottles each of 0%IS, 0.0032%IS, 0.01%IS, 0.1%IS, 1%IS, 10%IS x 4 mL

Xpert BCR-ABL 0%IS*

Catalog No.: C13112-5

Kit Contains: 5 bottles of 0%IS x 4 mL

Xpert BCR-ABL IS p210 Linearity Panel C207*

Catalog No.: C207

Kit Contains: 2 bottles each of 0.01%IS, 0.1%IS, 1%IS, 10%IS, 20%IS, 50%IS x 4 mL

**FDA-cleared IVD & traceable to WHO*

Xpert BCR-ABL1 p190 Control Panel C183

Catalog No.: C183

Kit Contains: 2 bottles each of 0%, 0.02%, 0.1%, 1%, 10% x 4 mL

Xpert BCR-ABL1 p190 0%

Catalog No.: C17612-5

Kit Contains: 5 bottles of 0% x 4 mL

COMING
SOON

Xpert BCR-ABL1 p190 Linearity Panel C208

Catalog No.: C208

Kit Contains: 2 bottles each of 0.008%, 0.02%, 0.1%, 1%, 10%, 20% x 4 mL



NPM1

Xpert NPM1 Control Panel C194v1.1

Catalog No.: C194v1.1

Kit Contains: 2 bottles each MUT 0%, MUT A 1%, MUT B 1%, MUT D 1%, and MUT A 5% x 4 mL

Xpert NPM1 MUT 0%

Catalog No.: C19512-5

Kit Contains: 5 bottles of 0% x 4 mL

COMING
SOON

Xpert NPM1 Linearity Panel C225

Catalog No.: C225

Kit Contains: 2 bottles each MUT A 0.1%, MUT A 1%, MUT A 20%, MUT A 100%, and MUT A 450% x 4 mL



PML-RARA

Xpert PML-RARA Control Panel C215

Catalog No.: C215

Kit Contains: 2 bottles each 0%, bcr1 1%, bcr2 1%, bcr3 1%, and bcr1 20% x 4 mL

Xpert PML-RARA 0%

Catalog No.: C21612-5

Kit Contains: 5 bottles of 0% x 4 mL

COMING
SOON

Xpert PML-RARA Linearity Panel C221

Catalog No.: C221

Kit Contains: 2 bottles each bcr1 0.1%, bcr1 1%, bcr1 20%, bcr1 100%, and bcr1 450% x 4 mL

