



Quality Controls for Respiratory Testing

Support confident clinical decisions by ensuring accurate results with **trusted molecular quality controls**. Our comprehensive respiratory control solutions help maintain accuracy, consistency, and compliance—test after test. MMQCI controls deliver proven, reliable performance **across a broad range of leading molecular testing platforms**.



Comprehensive & Multiplex

Covering a broad range of pathogens such as flu, RSV, SARS-CoV-2, adenovirus & more in only a few tubes!



Wide Platform Range

Ensure consistent, dependable results throughout a wide range of molecular testing platforms



Trusted worldwide

Proven performance recognized across the globe by leading healthcare professionals and labs



Convenient & Ready-to-Use

Available in liquid form, our controls can be used as provided, without reconstitution or pooling.



Reliable & Accurate

Proven lot-to-lot consistency protects accuracy and helps identify shifts or trends.



Innovative and Safe

Patented DNA/RNA controls contain no human or animal material, ensuring safe handling



Product Name**Catalog No.****Configuration****Storage****BIOFIRE**

BioFire® RP2.1/RP2.1plus Control Panel M441	M441	12 x 300 µL	Frozen
FilmArray® Pneumonia/Pneumoniaplus Control	M340	12 x 120 µL	Frozen
SPOTFIRE® RSP Pos & Neg Controls	M425	12 x 300 µL	Room Temp.
SPOTFIRE® RSP Positive Control	M42638	6 x 300 µL	Room Temp.
SPOTFIRE® RSP Negative Control	M42738	6 x 300 µL	Room Temp.

Roche

eplex RP Control M306	M306	10 x 200 µL	Frozen
eplex RP Negative	M31118-6	6 x 200 µL	Frozen
eplex RP2 Control M451	M451	6 x 200 µL	Refrigerated
<i>eplex RP3 Control M530</i>	<i>M530</i>	<i>12 x 200 µL</i>	<i>Refrigerated</i>

**QIAGEN**

QIAstat-Dx RP SARS-CoV-2 Control Panel	M360	12 x 300 µL	Frozen
QIAstat-Dx RP SARS-CoV-2 Negative	M36321-6	6 x 300 µL	Frozen

Diasorin

Verigene® RP Flex Control Panel M270	M270	12 x 200 µL	Frozen
<i>LIAISON PLEX RP Control Panel M518</i>	<i>M518</i>	<i>12 x 300µL</i>	<i>Refrigerated</i>

**Cepheid**

Xpert SARS-CoV-2 Control Panel M430	M430	6 x 300 µL	Frozen
-------------------------------------	------	------------	--------

Curetis

Unyvero LRT/Pneumonia Control Panel M404	M404	12 x 180 µL	Refrigerated
--	------	-------------	--------------

All other platforms and LDTs

Birlinn™ CoV2-Flu-RSV Control Panel M509	M509	12 x 800 µL	Refrigerated
--	------	-------------	--------------

