

NEW

Infectious Disease Testing on the Diasorin LIAISON PLEX®



MULTIPLEX

Monitors the amplification and detection of all 19 pathogens included in the LIAISON PLEX® Respiratory Flex Assay with only 2 tubes, saving costs of cartridges used



READY & SAFE TO USE

No reconstitution or pooling required. Contains carefully crafted genome segments that are 100% non-infectious and safe for the lab to handle.



TRUE NEGATIVE

Includes a true negative control that contains no target template, able to detect low-levels of contamination.



VERSATILE

Can be used for validation, verification, routine monitoring, and operator training, easy to include in your IQCP.



LIAISON PLEX RP Control Panel M518

Catalog No.: M518

Kit Contains: 12 tubes x 300 µL

4 each of Pos A, Pos B, and Negative controls

